

FIST-CROSS BELT SLING



HOW TO PLACE



Ref: FI-410
Ref: FI-410P



Introduce the hand of the injured arm in the cushioned wristband.
Put the hand at the desired height of the body and lean the band on the shoulder.



Turn the band round the back, across one's chest.

IN CASE OF ARM OR HAND INJURY

In case of arm or hand injury, it is recommended its fixation to the body in order to avoid involuntary hits.

Our product FI-410 is a cushioned wristband attached to a band made of different materials welded together by using an ultrasonic process to which it is possible to attach the velcro at the end.



Using this design the patient does not have to support the weight of the arm permanently with his neck.

With the forearm folded, add the band by using the velcro where it is supposed to.



Distributor

2mil artesanía sanitaria, s.l

✉: c/Amadeo Vives nº 8 - 28290 LAS MATAS (Madrid) SPAIN

☎ & FAX: +34-916301254

ER: +34-606957883

E-Mail: 2mil@2mil.com

www.2mil.com