

# NECK BRACES

CE




Ref: CO-050  
CO-100-I


Before wearing read the following considerations:

- The duration of the treatment should be decided by a doctor exclusively
- Expressly respect the daily duration of use of this product following your doctors indications
- Don't use any medications of "external use" in the area where applying this product without the authorisation of your doctor
- Use this product with the appropriate size.
- Recommended only wear during the day, removing at bedtime or when receiving physiotherapy treatments.

## CERVICAL NECK BRACES OF A MEDIUM DENSITY



Children's sizes – Large and Small (babies).



Adult's sizes – Large, Medium and Small.

This 2miil product is foam composed of medium-density 100% cotton lining and fixed by Velcro, allows mental and occipital support, so that limits the mobility of the cervical spine allowing a relaxation of the muscles. At the same time decreasing the pressure exerted by the weight of the head on the cervical vertebrae, which reduces the pressure at the disk.



Distributor

**2miil** artesanía sanitaria, s.l

✉: c/Amadeo Vives nº 8 - 28290 LAS MATAS (Madrid) SPAIN

☎ & FAX: +34-916309208 / 1254

ER: +34-606957883

E-Mail: 2miil@2miil.com

www.2miil.com