

# THORAX BELT



HOW TO TIGHTEN THE SLINGS



Ref: SU-225

TO BE USED ON LAYING DOWN PATIENTS.




The belt may be fastened to the bed by the straps provided.

The patient is easily positioned as the belt fastens from the back to the thorax.

Once closed, Velcro fastens it.

The fastening belt is attached in such a way that the patient can move from side to side.


1st INTRODUCTION:




2nd INTRODUCTION:



TIGHTEN:



FASTEN:



Distributor

**2mil** artesanía sanitaria, s.l

✉: c/Amadeo Vives nº 8 - 28290 LAS MATAS (Madrid) SPAIN

☎ & FAX: +34-916309208 / 1254

ER: +34-606957883

E-Mail: 2mil@2mil.com

www.2mil.com

