

TRUNK RESTRAINER



Ref: SU-220
SU-222



PROTECT PATIENTS FROM FALLS

Sedated or disoriented patients normally lose their sense of balance which, at times, means that they suffer from falls. To avoid this, the nursing staff usually restrain the patients using sheets, bandages or anything else which comes to hand which may be used for this purpose.

2mil's trunk retractor has been designed to restrict the patient movements to those permitted by the nursing staff for each individual case. There is available an Extra size, which is supplied with a longer belt to accommodate larger patients. (1)


The belt is made using an ultrasonic process and it therefore has no seams and skin irritation through rubbing is avoided.

It is extremely soft and can be placed directly onto the skin.

The restraining belts are highly resistant and they are sufficiently long for any kind of tying system to be used. Likewise, different levels of mobility may be permitted.

Your salesman will inform you as to the best way of using this equipment.

RESTRAINING OF PATIENTS IN BEDS OR CHAIRS



- (1) SU-220: Trunk Restrainer
SU-222: Trunk Restrainer Extra



2mil artesanía sanitaria, s.l

✉: c/Amadeo Vives nº 8 - 28290 LAS MATAS (Madrid) SPAIN
☎ & FAX: +34-916309208 / 1254 ER: +34-606957883
E-Mail: 2mil@2mil.com www.2mil.com

Distributor